



Outline

Vacuum desiccator provided by the manufacturer of the vacuum pump.

Standard/Size

Capacity: 3.3L

Inside dimension: Bottom I.D.165Xl.D.145XHeight175

Pressure: Guaranteed up to 13.3Pa

Base: 4 ports (Sealing plug X2, Hose port suction pipe (dia.9))

Cover material: Polycarbonate

Base material: Aluminum

O-ring material: FPM

Weight: 3.0kg

Overall dimension: 225(W)X213(L)X235(H)

Features and Cautions on use

Usable from low vacuum to medium vacuum.

Easy to see inside and clear visibility.

Can be used for multiple purposes with four connection ports.

Two valve sets are standard equipment.

There is lock function and is reliable.

Realizes high sealing performance and high strength with a metal base.

アルバック機工株式会社

〒881-0037

宮崎県西都市大字茶臼原291-7

TEL:(0983)42-1411(代)

FAX:(0983)42-1422